

LG Quad Combo™ Quick Start Guide

Step 1: Take out all contents in the box and Insert the Battery



Step 2: Connect both lead wires to the top of the unit



Step 3: Connect the lead wires to the electrodes, place all 4 electrodes on effected area



Step 4: Turn the Device “On” by pressing the silver round bottom on the front of the device



Step 5: Press “M” to Select IF (interferential), EMS (E-Stim), Tens or Micro (Microcurrent). For which waveform and program to choose, please consult your prescribing healthcare provider.



Step 6: Select the up arrow on the left until the electrodes connected to that lead wire are strong but comfortable. If using all 4 electrodes, repeat this step using the up arrow on the right side.



Step 7: Enjoy the pain relief and muscle strengthening settings available with this product. For more specific information, please consult the user manual or your Healthcare Provider.